

User Manual


GEW Series

Digital Weighing Scale

Cautions

- Upon receiving the scale, please take at least 8 hours to recharge the battery before you start to operate this scale.
- Please recharge the battery after a long period of storage (more than one month). It takes 8~10 hours to fully recharge the battery.
- Please assemble the platter before powering on the scale, especially for the scale with 1.2kg / 3kg / 6kg max capacity.
- Do not keep the weight placed on the platter for a long period.
- Keep the scale away from the environment with high temperature, high humidity, heavy press, and heavy bump.
- Always make sure the scale is located in a flat and plane surface.

Dimensions



Symbols on LCD

■ Zero Symbol (O) :

O appears on the left-bottom corner of LCD, and indicates the scale is at zero-point.

■ STABLE Indication (▼) :

▼ cursor appears on the bottom edge of WEIGHT LCD (above of **STABLE** printed on the display-overlay), and indicates the scale is at stable situation.

■ Indication of Weighing Range (▼) :


▼ cursor appears on the bottom edge of WEIGHT LCD (above of **W1** or **W2** printed on the display-overlay), and indicates the present weighing range to be **Weighing Range 1 (W1)** or **Weighing Range 2 (W2)**.

!! Notice: *The indication of weighing range is only available for the scale with multi-range (2 ranges).*

■ NET Weight Symbol (Net) :


Net appears on the bottom edge of LCD, and indicates the tare weight has been deducted.

■ Battery Symbol () :



 appears on the bottom edge of LCD, and indicates the battery is running in low-voltage status. This is for reminding the user to recharge the battery.

Operation Instructions

■ Backlight (optional) :

Press and hold →  ← key till the scale sounds **beep** to enable / disable the Backlight.

■ Status of the Battery :

When LCD indicates  and the scale sounds **beep beep beep** every thirty-minute, it means the battery is running in low-voltage status. LCD will indicate  continuously for reminding the user to recharge the battery. If the battery is still running in low-voltage status without being recharged, LCD will indicate **-OFF-** and the scale will sound **beep** repeatedly. Please power off the scale and recharge the battery immediately.

■ UNIT:

Select the weight mode: **kg** / **pcs**.

■ →0← :

Pressing this key to make the scale be at zero point. LCD will indicate **0** on the left-bottom corner.

■ ←T→ :


Pressing this key to deduct the tare weight. Put the package or container on the platter and then press [**TARE**] key. LCD will show **0** and **Net**.

■ Under **pcs** (simple counting) mode:

!! The suggested unit weight of the measured subject should be heavier than **the division of this scale**.

1. Put 100 pieces of the measured subject (with the same unit weight) on the platter.
2. Press and hold [**UNIT**] key till LCD shows **S= 10**, and then quickly (within 3 seconds) press [**UNIT**] key again to select **S=100** (10 · 20 · 50 · 100 could be selected for the sample-quantity of the measured subject).
3. Wait till LCD shows **100**. The scale has finish calculating the averaging unit weight of the measured subject. Now the simple counting function is ready. This scale could be used as a simple counting scale.

RS232 Interface Connection: *<optional device>*



RS232 Wiring Configuration:

Baurate: **9600**

Parity: **None**

Stop Bit: **1**

Data Bit: **8**

Recharge the Battery

Power off the scale and connect the power cord to an AC outlet. CHARGE LED will indicate the ongoing status of the battery. It takes about **8~10 hours** to fully recharge the battery.

CHARGE LED:

Color of LED	Status of the battery
RED	Initial Connection
ORANGE	Charging
GREEN	Fully-charged

Power Requirements

- AC: 220V, 50/60 Hz
- DC: 6V/25mA Rechargeable Battery;
P=0.2W(max)

We do more than you expect
